

aeroGAIN-BASE-I.I

High power ytterbium fiber gain module

- Truly single mode polarization maintaining system
- Excellent pointing stability
- Robust industrial construction
- 10 μm or 15 μm step-index fiber input
- Easy thermal management
- Long lifetime



Applications

Ultrafast fiber lasers




The *aero***GAIN**-BASE-1.1 is a high performance ytterbium fiber gain module designed for industrial manufacturers of pulsed fiber lasers, and is also suited as an easy entry into high power scientific setups. The module is equipped with a 10 or 15 μ m step-index standard single-mode fiber input that can easily be spliced to a seed source. The gain medium is our world-renowned DC-200/40-PZ-Yb fiber providing the largest single mode MFD in the industry. The gain fiber is heatsunk to the aluminum base-plate that can be clamped to e.g. a water chilled plate or an air-cooled heat sink. The output end of the module is equipped with a large AR coated endcap that provides mode expansion and reduces reflections. Excess pump light is removed by the integrated residual pump dump.

For optimal performance, the *aeroGAIN*-BASE-1.1 is designed for counter propagating pumping through the output endcap. The end-cap and the last part of the gain fiber is mounted in a water cooled housing to ensure maximum performance and lifetime.

125/10 step-index input seed fiber	Residual pump dump	40 μm core single-mode gain fiber	Water cooled section	End-cap
Model	YI	Yb Fiber Length		ignal
aeroGAIN-BASE-1.1		3.0 m	1064 nm	

All modules are assembled and tested in cleanrooms, and come mounted on a tooling plate to ensure the endcap is clean. This tooling plate can also can be used as a mount in the laboratory.

The aeroGAIN-BASE design has proven its industrial 24/7 reliability through a significant number of long term tests including continuous 25,000 hours operation at 55 W output power. Lifetime in specific OEM systems depends on pump power, pump wavelength, cooling efficiency etc. and NKT Photonics is available for supporting system design.





Other aeroGAIN products

aeroGAIN-ROD

When a 40 μm core is not big enough our aeroGAIN-ROD modules deliver the ultimate gain solution. With MFDs in excess of 65 μm , ROD systems sits at the very top of the amplifier chain reaching power levels normally only found in DPSS systems but with the benefits of a fiber waveguide and the efficient ytterbium material system.



NKT Photonics A/S (Headquarters)

Blokken 84, 3460 Birkerød, Denmark

Phone: +45 4348 3900 Fax: +45 4348 3901

NKT Photonics GmbH

Schanzenstrasse 39, Bldg D9-D13 51063 Cologne, Germany Phone: +49 221 99511-0

Phone: +49 221 99511-0 Fax: +49 221 99511-650

NKT Photonics Inc.

Office 23, 4400 Route 9 South Freehold, NJ 07728, USA Phone: +1 732 972 9937 Fax: +1 732 414 4094



Specifications

Optical

Seed input	
Signal wavelength	1060 – 1070 nm
Recommended signal input power	> 500 mW for high power operation
Recommended pulse duration	fs, ps or ns pulses shorter than 2 ns
Signal input fiber	- 10 μm core PM 125 μm / 250 μm or - 15 μm core PM 250 μm / 350 μm
Signal Output	
Max signal gain	< 20 dB
Rated output power	75 W
M ²	≤1.3
Mode field diameter	31 ± 2 μm
PER	≥ 15 dB
Typical optical efficiency*	> 70 %
Typical core to clad power ratio*	> 96 %
Pump Input	
Pump center wavelength	976 ± 2 nm
Maximum pump power P _p	100 W @ fiber facet
Recommended pump type	Fiber delivered 200/0.22 (max NA<0.55)
Pump cladding diameter	200 ± 2 μm

 $[\]mbox{*}$ Evaluated with 2 W input power at 1064 nm and 75 W output power.

Mechanical

Base plate dimensions (HxBxW)*	See drawing
Weight	1.8 kg (without tooling plate)
Length of input pigtail	1 m
Output end facet angle	0 degree
Endcap length / diameter	6 mm / 7 mm** - AR coated

^{*} The system is shipped on a larger tooling plate that can also be used for mounting the module during test

Water cooling

Cooling flow for base	4 ± 2 liter/minute
Cooling flow in the tube	0.2 - 0.5 liter/minute
Cooling water temperature	25 ± 5°C
Max temperature of base-plate	35 °C

All NKT Photonics products are produced under our quality management system certified in accordance with the ISO 9001:2015 standard.







^{**} Open aperture